

# Drying Cabinets

After cleaning and disinfecting medical **Hospital CSSD Equipment**

devices, dry them in your CSSD LuxMed® drying cabinet: Single-door or pass-through, 4+4 shelves, accessories for various needs.

## Application:

Hot air sterilizer is used to sterilize small surgical instrument, glass, petri dishes etc.

This device is suitable for medical and veterinary clinics, hospitals, pharmacies, health care centers, and laboratories.

## Features:

- Adjustable sterilization time.
- Over-temperature alarm.
- Top air outlet with adjustable air volume.
- Heating up fast, forced-air convection, fast sterilization.
- High precision microprocessor temperature control; 304 stainless steel inner chamber.
- Fan inside, fully sealed inside, no air circulation with external air, temperature rise fast.



## Models:

Model	LMS-T25	LMS-T50	LMS-T70	LMS-T50	LMS-T50
Capacity	25L	50L	70L	50L	50L
Temp. Range	RT+10~250°C				
Temp. Precision	0.1°C				
Temp. Fluctuation	±0.5°C				
Ambient Temp.	5~40°C				
Timing Range	1~9999min				
Shelves No.	2 pcs				
Consumption	500W	750W	1050W	1500W	2000W
Power Supply	AC110/220V±10%,50/60Hz				
Internal Size(W*D*H) mm	300*300*270	420*350*350	4400*375*470	550*350*520	600*550*570
External Size(W*D*H) mm	590*470*440	710*520*520	540*550*870	840*540*720	890*730*770
Packing Size (W*D*H) mm	710*600*550	810*640*600	670*660*1040	950*680*880	1020*890*930
Gross Weight(kg)	45	50	75	80	100

M: +91 8130383561, 9350831213 | E: cto@luxmed.in

SCAN & SEND MESSAGE :



SCAN & FOLLOW US ON :

