

Drying Cabinets HEPA-filtered

Hospital CSSD Equipment

LuxMed[®]
Healthcare Redefined

LuxMed® presents state-of-the-art

means for drying by HEPA-filtered air

Special attention has been paid to the flexibility of use of the large volume chamber by providing purpose designed shelves and inserts for accommodating a variety of medical devices.

The Programmable Logic Control system with advance User Interface & Display panel allows fixed and tailored cycles and easy follow-up of ongoing process status.

Single Door or Pass-through, left/right hinged door(s)



Features:

- Sound & heat insulated
- Heat exchanger
- Single door or pass-through, left/right hinged door(s)
- Adjustment of the interval between retractable shelves
- Advance control panel & info display
- Temperature adjustment range 30-80 C,
- Time adjustment range 20 -180 minutes / non-stop
- Cooling time range 1-30 minutes
- Temperature & remaining time indication
- HEPA -filter condition indication
- Double overheating protection

Models:

Model	LMS-T25	LMS-T50	LMS-T70	LMS-T50	LMS-T50
Capacity	25L	50L	70L	50L	50L
Temp. Range	RT + 10~250°C				
Temp. Precision	0.1°C				
Temp. Fluctuation	±0.5°C				
Ambient Temp.	5~40°C				
Timing Range	1~9999min				
Shelves No.	2 pcs				
Consumption	500W	750W	1050W	1500W	2000W
Power Supply	AC110/220V ± 10%, 50/60Hz				
Internal Size(W*D*H) mm	300*300*270	420*350*350	4400*375*470	550*350*520	600*550*570
External Size(W*D*H) mm	590*470*440	710*520*520	540*550*870	840*540*720	890*730*770
Packing Size (W*D*H) mm	710*600*550	810*640*600	670*660*1040	950*680*880	1020*890*930
Gross Weight(kg)	45	50	75	80	100

M: +91 8130383561, 9350831213 | E: cto@luxmed.in

SCAN & SEND MESSAGE :



SCAN & FOLLOW US ON :

